




Product List WESTERNBLOT and EUROLINE-WB



Type	Antigens
Liver with SLA/LP (EUROLINE-WB)	Liver extract plus rec. SLA/LP
HEp-2 (EUROLINE-WB)	Full extract plus native SS-A and Rec. Ro-52
Mi-2	Cell extract
Neuronal AB (EUROLINE-WB)	Cerebellum extract plus rec. Hu and rec Yo
Borrelia with VisE (EUROLINE-WB)	Full extract plus rec. VisE
Borrelia burgdorferi sensu stricto	Full extract
Borrelia garinii	Full extract
Borrelia afzelii	Full extract
Y. enterocolitica	Full extract
HSV 1 und 2 (EUROLINE-WB)	HSV 1 full extract plus native purified HSV 2 gG2
H. pylori	Full extract
EBV	Full extract
Treponema pallidum	Full extract
Echinococcus granulosus	Full extract
CMV	Full extract

version (07/03)